

## **Uni-Grip™ Universal Disposable Radiography Holder**

Taking x-rays, whether with film, phosphor plates or digital sensors, is simpler with the Uni-Grip Universal Disposable Radiography Holder. Unlike other disposable holders on the market, which require multiple units to perform a full mouth x-ray series, you only need one Uni-Grip for anterior, posterior and bitewings radiographs. Accommodates any #1 or #2 size digital sensor, film or phosphor plate. Uni-Grip is fast, accurate and its compact size and rounded corners are easy on the patient.

Please visit these links for additional information.

[http://www.rinncorp.com/catalog\\_sensor\\_unigrip.shtml](http://www.rinncorp.com/catalog_sensor_unigrip.shtml)

<http://www.rinncorp.com/pdf/UniGrip%20Flyer.pdf>

For training guide go to this link <http://www.unigrip.com/>



Item No. 55-0050 - Uni-Grip™ Disposable Holder for Digital Sensors (box-50 Uni-Grips + 50 Digital Sensor Barriers)



Item No. 55-0250 - Uni-Grip™ Combo Kit for Digital Sensors (2 55-0050 kits + 1 55-0100 Uni-Grip™AR)

Item No. 55-0100 - Uni-GripAR Positioning Arm and Ring for Uni-Grip Universal Disposable Radiography Holders