

Suni Windows Compatibility Matrix

04-Oct-11 mpr

Product	Interface	Windows XP (32-bit)	Windows XP (64-bit)	Windows Vista (32-bit)	Windows Vista (64-bit)	Windows 7 (32-bit)	Windows 7 (64-bit)	
SuniRay (CMOS sensor)	Integrated USB	Yes	Yes	Yes*	Yes*	Yes	Yes	
Dr. Suni Plus (CCD sensor)	USB 2000	Yes	Yes	Yes*	Yes*	Yes	Yes	
Dr. Suni Plus (CCD sensor)	USB 1200 / 1600 /1800	Yes	No	Yes	No	No	No	
Legacy/metal (CCD sensor)	USB 1200 / 1600 /1800	Yes	No	Yes	No	No	No	
SuniCam	Integrated USB	Yes	No	v1.8 drivers required		No	v1.8 drivers required	No
SuniCam II	USB Docking Station	Yes	Yes	Yes	Yes	Yes	Yes	
Suni 3D (all models)	Proprietary CAT5	Yes	No	Yes	No	Yes	No	

* Install using Windows XP instructions