

XCP-ORA™

ONE RING & ARM POSITIONING SYSTEM

For **Film** and **Sensors**

Blue and yellow pins for anterior and posteriors



Red pins for bitewings



Yellow: Posteriors




Red: Bitewings




Blue: Anteriors



Anteriors



Posteriors



Horizontal & Vertical Bitewings

The new XCP-ORA™ One Ring & Arm positioning system reduces the number of components needed for accurately positioning radiographic media with RINN XCP® or XCP-DS™ biteblocks for Anterior, Posterior or Bitewings to only two pieces.

- **Autoclavable** for low cost-per-use
- **Lightweight construction** helps patients maintain positioning
- **Patent-pending Ring & Arm Design** simplifies your positioning procedure with less parts
- **Patented color coding** for faster, easier, more accurate positioning
- **More economical** than separate rings and arms

550771 XCP-ORA™ One Ring & Arm, 1-pack MSRP \$104.99

550772 XCP-ORA™ One Ring & Arm, 2-pack MSRP \$189.99

XCP® and Rinn® are registered trademarks, and XCP-DS™ and XCP-ORA™ are trademarks of DENTSPLY International, Inc.

Patent pending.

©2009 DENTSPLY Rinn MK201 3.09

DENTSPLY
RINN

A Division of DENTSPLY International, Inc.

1212 Abbott Drive

Elgin, Illinois 60123-1819 USA

800 323 0970

847 742 1115

847 742 1122 Fax

www.rinncorp.com

www.unigrip.com