

USING RINN XCP or ORA holder with your SuniRay sensor

Suniray RINN Bite blocks part numbers are:

BiteBlocks	Part number
#1H post	552111
#2H Post	552112
#1V Ant	552113
#2V ANt	552114
#1V Endo	552119
#2H Endo	552120
#1H Bitewing	550623
#2H Bitewing	550624
#1V Bitewing	550625
#2V Bitewing	550626

You can use the RINN ORA RING and ARM with these bite blocks. Visit this link for additional information.

http://www.rinncorp.com/pdf/MK201_XCPORA.pdf

Or you can use RINN XCP-DS Rings and ARM with these bite blocks. Visit this link for additional information.

http://www.rinncorp.com/pdf/DentsplyRinn_2009Catalog.pdf

Visit this link for training on these sensor holders.

<http://www.rinncorp.com/modules.shtml>